



 Convaid

# Convaid Cruiser

Lightweight, Compact-Folding, Convenient Fixed-Tilt Buggy



# Convaid Cruiser

The Convaid Cruiser is a popular 30° fixed-tilt pediatric buggy that is convenient, lightweight, and compact folding. The Cruiser is made to order to meet each individual requirements. It comes with Convaid's proprietary self-tensioning seating system, Self-Tensio®, which provides increased seating tolerance and can help promote neutral pelvic positioning. It provides growth capability through its seat depth extension system.



# Lightweight, Compact-Folding, and Convenient

The Convaid Cruiser is a 30° fixed-tilt buggy that supports postural control through upper body positioning.



## Tension adjustable back upholstery

Adjustable straps allows individual adjustment.



## One-piece height adjustable push handle

Two-piece push handle option (standard for size 18)



## Self-Tensio®

Self-Tensio® is Convaid's unique self-tensioning seating system which provides pelvic positioning similar to a solid seating system in an umbrella folding rehab buggy.



## Tool-less upholstery removal

No tools required to remove for cleaning

## Standard

- 30° Fixed Tilt
- One Piece Height Adjustable Push Handle (CX 10/12/14/16)
- Two Piece Push Handles (CX 18)
- Tension Adjustable Back
- Tool-less Upholstery Removal
- Tool-less Adjustable Seat Depth Growth of 13mm
- Zippered Pockets in Seat & Back Upholstery (Cordura Only: CX10/12/14/16)
- 2-Point Positioning Belt (Not available in EU)

- Self-Tensio®
- Elevating Legrest
- Removable, Height Adjustable Swing-Away Footplates
- 7,5 x 2" Front, 12,5 x 2" Rear Solid Tire (CX10/12/14/16)
- Quick Release Rolko Rear Wheels (CX10/12/14/16)
- Foot Operated Wheel Lock (CX10/12/14/16)
- Toggle Wheel Locks (CX18)
- Silver Powder Coated Frame
- Made in the USA

## Options

### Sizes

- 10
- 12
- 14
- 16
- 18

### Seat Depth Tubes

- 51mm seat depth tube extensions
- 127mm seat depth tube extensions
- Custom seat depth extension 76mm

### Upholstery

- Cordura
- Heavy Duty Cordura

### Tires

- 12,5" Rear Solid Knobby Tire
- 7,5" x 2" Front, 12,5" Rear Air Knobby Tire
- 8" x 2" Front Ribbed, 12,5" Rear Solid Knobby Tire
- 8" x 2" Front Ribbed, 12,5" Rear Air Knobby Tire

### Cordura Upholstery Colors



Princess Pink



Cherry Red



Apple Green



Forest Green



Aqua Blue



Royal Blue



Sassy Purple



Chocolate Brown



Panther Black

(limited range of colour options in specific countries)

### Frame Color



Cosmic Silver

### Frame Options



Scout



Size 18



## Accessories



**Extended Headrest Cover**  
- w/and without windows  
(Canopy)



**Headrest Cover (Canopy)**



**Rain Cover**



**Mosquito Net**



**Headrest Extension -  
Rectangular - Textilene**

(Not available in EU)



**Headrest Extension -  
Curved - Cordura**



**Occi-Headwings**



**Padded Headwings**



**Seat Reducer Insert**



**Full Torso Support Vest**



**Soft Adjustable Lateral  
Support Single Flap w/  
Scolio Strap**



**Soft Adjustable Lateral  
Support Double Flap w/  
Scolio Strap**

## Accessories



**H-Harness with padded covers**

3-Point Positioning Belt Required



**General Use Cushion<sup>\*1</sup>**

(Not available in EU)



**Align Cushion - Lateral Pelvis/Lateral Thigh Support<sup>\*1</sup>**

(Not available in EU)



**Position Cushion - Medial Thigh Support/Anti Thrust<sup>\*1</sup>**

(Not available in EU)



**3-Point Positioning Belt**



**Lateral Thigh Support - Abductor**



**Medial Thigh Support - Abductor**



**Flip-Up Armrests**



**Anatomic Back**

Heavy Duty Reinforced Upholstery Required



**Attendant Hand Brake**



**Upper Extremity Support Surface (Clear Tray)**



**Under Seat Storage Basket - Fixed Tilt**

<sup>^</sup> Crash tested Foam Tray with vinyl cover. Flip-Up Armrests required"

<sup>\*</sup>Manufacturer Installation Required

<sup>‡</sup>Incontinence Cushion Liner available

## Accessories




**Saddle Bags**



**Utility Bag**




**Travel Bag**



**Calf Panel**



**Angle Adjustable  
Footplates**



**Footplate Securement  
Strap**




**Castor Lock for 2" Wide  
Front Tires**

(Standard in EU)



**5-Point Harness with  
Padded Covers**

Crash Tested, WC19 Compliant  
(Cruiser 10+12 only)\*  
(Not available in EU)





















**Transit Anchors**

\*Manufacturer Installation Required

Limited range of accessory options in specific countries

# Measurements

Cruiser measurements		CX 10	CX 12	CX 14	CX 16	CX 18
	Seat Width (mm)*	254	305	356	406	457
	Seat Depth (mm)	152-279	203-330	254-381	279-406	356-533
	Seat Back Height (mm)	533	559	635	686	635
	Seat to floor (mm)	584	635	533	584	660
	Seat to Back Angle	85°/90°/95°	85°/90°/95°	85°/90°/95°	85°/90°/95°	85°/90°/95°
	Seat Angle	30°	30°	30°	30°	30°
	Overall Height (mm)	940	965	940	1016	965
	Overall Length (mm)	1219	1118	1245	1321	1397
	Overall Width (mm)	559	597	648	699	724
	Folded Height (mm)	432	432	432	432	432
	Folded Width (mm)	381	381	406	406	381
	Folded Length (mm)	1194	1219	1194	1346	1270
	Headrest Extension (mm)	152	203	203	203	254
	Shoulder Strap Height (mm)	279-483	305-508	381-584	432-635	406-559
	Seat to Footplate (mm)*	127-533	127-635	102-559	102-559	203-584
	Push Handle Adjustments (mm) °	711-1219	737-1168	711-1143	787-1219	N/A
	Weight Capacity Regular/Transit	34/30 kg	34/30 kg	45/45 kg	77/77 kg	113/77 kg
	Weight of Chair without Front Rigging	12 kg	12 kg	13 kg	15 kg	14 kg

\* Plus or minus 13mm due to fabric  
° Adjustable push handle is not available on CX 18 model

Convaio-Cruiser-BR1901-EUR-REV.00 2022-07-11 | Reserving our rights for printer's error discontinued products. Please find an updated edition on our website © 2022 Etac All rights reserved

For the latest news and continuously updated product information visit: [www.etac.com](http://www.etac.com)