



 **Convaid**

# Convaid Rodeo

Lightweight and Convenient Tilt-in-Space Buggy



# Convaid Rodeo

The Convaid Rodeo is a popular pediatric tilt-in-space buggy that is convenient, lightweight and folds easily. The Rodeo provides up to 45° of adjustable tilt that offers pressure relief and upper body support. Made in order to meet user's specific positioning needs and provides 5 years of growth capability through its seat depth feature. With 2 frame options available (with and without seating), the Rodeo can accommodate aftermarket and extensive seating needs.



# Lightweight, Compact Folding, and Convenient



## Up to 45° adjustable tilt

Provides pressure relief and gives the child a possibility to rest.



## Easily fold

Folding mechanism allows the whole chair, including seating, to fold into a tight (<457mm H) and lightweight (13-18kg) package which fits in the trunk of almost every car.



## Tension adjustable back upholstery

Adjustable straps allows individual adjustment.



## Zippers in seat for foam inserts

Allows for additional seat padding and support.

## Standard

- 5° to 45° of Adjustable Tilt
- One Piece Height Adjustable Push Handle
- Seat Depth Growth of 127mm
- Tool-less Upholstery Removal
- Zippered Pockets in Seat Bottom Upholstery
- 2-Point Positioning Belt (Not available in EU)
- Elevating Legrest
- One Piece Angle Adjustable Footplate
- Option for solid or air tires
- Quick Release Rolko Rear Tires
- Single Foot Activated Wheel Lock Bar
- Silver Powder Coated Frame
- Made in the USA

## Options

### Sizes

- 10 (Not available in EU)
- 12
- 14
- 16

### Transit Options

- 4 Black Transit Anchors
- Upper Extremity Support Surface Foam Tray


### Seat Depth Range Options

- Position 1: 178 - 305
- Position 2: 229 - 357
- Position 3: 279 - 405
- Position 4: 330 - 457

### Tire Options

- 7,5 X 2" Front Air Tires
- 12,5" Rear Solid Knobby Tires
- Castor Lock (EU Standard)

### Upholstery Colours



Panther Black



Cherry Red



Sassy Purple



Princess Pink



Royal Blue



Forest Green



Chocolate Brown



Aqua Blue



Apple Green

(limited range of colour options in specific countries)

### Frame Colour



Cosmic Silver

### Frame Options



Without seating (Not available in EU)



With seating and accessories



## Accessories



**Extended Headrest Cover**  
- with or without windows  
(Canopy)



**Headrest Cover (Canopy)**



**Rain Cover**



**Mosquito Net**




**Headrest Extension**



**Occi-Headwings**



**Padded Headwings**



**H-Harness**  
(Height adjustment harness by  
increment of 2")  
3-Point Positioning Belt Required



**Seat Reducer**



**Adjustable Firm Lateral  
Trunk Support -  
Contoured**



**Full Torso Vest**



**3-Point Positioning Belt**



## Accessories



**Attendant Hand Brakes**



**General Use Cushion<sup>\*1</sup>**

(Not available in EU)



**Align Cushion - Lateral Pelvis/Lateral Thigh Support<sup>\*1</sup>**

(Not available in EU)



**Position Cushion - Medial Thigh Support/Anti Thrust<sup>\*1</sup>**

(Not available in EU)



**Upper Extremity Support Surface (Clear Tray)**



**Upper Extremity Support Surface (Foam Tray)**

Crash tested Foam Tray with vinyl cover. Flip-Up Armrests required



**Medial Thigh Support (Abductor)**



**Lateral Thigh Support (Adductor) (not shown)**




**Solid Side Panels**



**Mesh Side Panels**



**Flip-up Armrests**



**Recline Option (90° to 110°)**

<sup>^</sup> Crash tested Foam Tray with vinyl cover. Flip-Up Armrests required"

<sup>\*</sup>Manufacturer Installation Required

<sup>1</sup>Incontinence Cushion Liner available



## Accessories



**IV Pole (height adjustable and collapsible)**



**Oxygen Tank Bag**



**Medical Necessity Storage Basket**


Under Seat Storage Basket - Tilt-In-Space (not shown)




**Utility Bag**



**Padded Elevating Legrest Cover**



**Elevating Legrest Extension - 51mm**



**A. Calf Panel  
B. Foot Positioners**



**Individual Angle Adjustable Footplates**



**Full Padded Footbox**  
One-piece footplate required



**Rear Anti-Tip Tubes (Pair)**



**Curb Tipper**  
(Standard in EU)














**Transit Anchors**

Limited range of accessory options in specific countries



# Measurements

Rodeo measurements		RD 10	RD 12	RD 14	RD 16
	Seat width (mm)	255	305	355	405
	Seat depth (mm)	180-305	230-360	280-4100	330-460
	Seat back height (mm)	530	560	610	635
	Seat to floor (mm)	510	530	560	560
	Recline adjustment (mm)	90° or 110°	90° or 110°	90° or 110°	90° or 110°
	Tilt adjustment (mm)	5-45°	5-45°	5-45°	5-45°
	Overall height (mm)	1090	1090	1120	1140
	Overall length (mm)	775	890	890	915
	Overall width (mm)	530	585	635	685
	Folded length (mm)	850	860	860	1015
	Folded width (mm)	530	585	635	685
	Folded height (mm)	445	430	460	460
	Headrest extension (mm)	150	200	200	200
	Shoulder strap height (mm)	330-480	380-530	410-560	460-610
	Footplate adjustment angle one piece angle adjustable footplate	±20°	±20°	±20°	±20°
	Seat to footplate one piece angle adjustable footplate	130-285	130-285	180-340	235-390
	Seat to footplate individual angle adjustable footplates	115-235	115-235	110-280	185-340
	Weight capacity regular/transit (kg)	34/30	34/30	45/45	77/59
	Weight of chair without front rigging (kg)	13	15	17	18
	Push handle adjustment (mm)	810-1230	880-1290	890-1320	930-1360

\*Reduced height when footplates are flipped

For the latest news and continuously updated  
product information visit: [www.etac.com](http://www.etac.com)